




## Wrist Wrap Application



Lay the Wrist Wrap flat on a table with the blue side up. Place your forearm on top of the Wrap with the palm facing up. Slip your thumb through the elastic thumb opening. Lay your hand palm-down if using on opposite hand.



Roll your hand over as shown. Securely fasten the hook strip around your forearm, then attach the hook strap around your hand. Be sure your thumb is outside of the Wrap.




Be careful not to over tighten the small hook strap around your fingers, so that when the Wrap inflates, it will not restrict blood flow.

**Note:** You may prefer to connect the hose before putting on the Wrap. For maximum cooling, place Wrap against bare skin.



CoolSystems, Inc.  
1201 Marina Village Pkwy Suite 200  
Alameda, CA 94501  
1.888.GAMEREDY

EC REP  
EMERGO EUROPE  
Molenstraat 15  
2513 BH, The Hague  
The Netherlands



Thank you for choosing the Game Ready System. To assure you have the best possible experience, please be sure to read the entire User's Manual prior to the use of this product.

### Warnings

- Improper placement or prolonged use of Game Ready could result in tissue damage such as frostbite.
- During the course of Game Ready therapy, patients should monitor the skin surrounding the treated region or the digits of the extremities of the treated limb for any burning, itching, increased swelling, or pain. If any of these signs present, or any changes in skin appearance occur (such as blisters, increased redness, discoloration, or other noticeable skin changes), patients are advised to discontinue use and consult a physician.
- The Game Ready Wrap is not sterile; do not place directly against open wounds, sores, rashes, infections, or stitches. May be applied over clothing or dressing.
- The Game Ready wraps are available in multiple configurations but are not intended for all possible orthopedic uses. Use caution when selecting a particular wrap for use on anatomic regions other than that for which the wrap is indicated and closely supervise the patient's status. For example, the Ankle Wrap is not designed for use on the toes and the Back Wrap is not designed for use in the abdominal region.
- If pressure is uncomfortable, pause unit, and restart with lower pressure.

### Care:

Hand wash the Wrap in warm water with a mild detergent or anti-bacterial soap. **Do not machine wash.**

### Warranty Information:

In case of manufacturer defect, Wrap may be returned within 1 year from date of purchase. See warranty card included with your Wrist Wrap.